



Ultra low temperature freezers -40°C and -86°C NEXUS Series





# OUR EXPERIENCE in biomedical cold storage



For over 85 years Angelantoni has been producing a full range of freezers for the long-term preservation of biological samples used in scientific laboratories, in pharmaceutical research, in hospitals and in any other application where it is necessary to preserve material of human, animal, plant origin or microorganisms at low / ultra-low temperatures.

Angelantoni's experience in the refrigeration market, combined with the continuous search for solutions which can provide its customers quality products, has led to the development of NEXUS, a new state of the art line of ultra-freezers with functions developed and tested by biomedical professionals.

The innovative handle and the new generation controller contribute to the **NEXUS** line's top level performance.





# Ultra low temperature freezers -40°C and -86°C **NEXUS** Series









# **NEXUS:** Product Features



#### **Maximum protection**

The AISI 304 Stainless steel internal chamber, together with the sophisticated external powder coated epoxy paint process, ensures maximum corrosion protection and long life.



#### **Reduced ice formation**

Inner doors with a variable configuration (from 3 to 5) and new concepts of pressure equalization valve reduce ice formation. The door's insulation efficiency is also guaranteed by a system of flat hinges which improve the seal of the gaskets.



#### **Outside communication**

Standard 23mm access port for the insertion of additional probes.



#### Interior space optimization

All ULT freezers are equipped with a system of height-adjustable shelves (drawers optional) allowing for the maximum utilization of the interior space, as well as the appropriate placement of the storage racks and cryogenic boxes.



#### Easy opening system

The main door opening and closing is characterized by a new, effective, ergonomic handle that allows for single handed operation. Standard key lock or optional BIOGUARD<sup>®</sup> lock system.



#### **Environmentally friendly**

At ANGELANTONI we consider ourselves to be "eco-friendly", so we use HCFC-free gas in all our products.









# J

#### Easy cleaning

The washable and reusable air filter is located in the lower front face of the ULT freezer, allowing easy access. Rounded corners of the internal chamber simplify cleaning.

# 0

#### Low noise levels

The cabinet and cooling system are designed to minimize noise, so the ULT freezers can be placed in laboratories or hospitals without disturbing operators.

## 

#### More safety in the accesses

The security of samples can be increased with the optional BIOGUARD<sup>®</sup> system, which provides complete access control.

With this option, the door is equipped with an electric lock that can only be opened using authorized cards. This ensures total traceability of accesses.

#### Easy handling

All vertical and horizontal ULT freezers are equipped with pivoting castors, making it easy to move them when necessary. Adjustable feet are provided for easy levelling.



#### **Freezing uniformity**

The serpentine cooling system ensures and constant cooling, guaranteeing the internal temperature stability and uniformity.

$\sim$	
<b>÷</b>	

#### **Maximum connectivity**

The electronic management of all ULT freezers is provided by the latest generation controller that optimizes performance and functionality and allows for simple remote adjustment and monitoring of the freezer.







Ultra low temperature freezers -40°C and -86°C **NEXUS** Series



# **NEW** ergonomic handle

### **EASY ACCESS TO YOUR SAMPLES**



Following an extensive study and design process, **NEXUS** ULT freezers have been equipped with a **new door handle** incorporating modern design and ergonomics. The **lever-operated** opening and closing mechanism, reduces the effort during opening and generates maximum pressure with minimum effort when closing. This makes it possible to increase the effectiveness of the insulation from the outside, reducing ice formation. The front position of the handle allows the ULT freezer to be positioned alongside other equipment, maximizing space usage.

In order to improve the mechanical performance, the handle also has a servo assisted system that totally eliminates friction, allowing the opening and closing to be done **with just one hand.** 

## More practical, improved performance.

In the **NEXUS** ULT freezers the handle is easily gripped and functional, providing the perfect balance between the Made in Italy design and ergonomics.







# LATEST GENERATION controller

### A LEADING COMPONENT FOR ITS RELIABILITY, PERFORMANCE AND EASY USE



**NEXUS** ULT freezers are fitted with the latest generation **ACP7** controller, developed after years of experience and utilizing the latest technologies. It is capable of guaranteeing high performance, maximum safety and ease of use.

The **controller** has a **user friendly** interface which allows the monitoring of functions and change of settings.

Access to parameters is password protected with **3 access levels**: user, service and administrator.

The sophisticated **display uses advanced technology** and is easy to use, thanks to **intuitive icon surfing**, the **touch screen** and the **multilingual instructions.** Its large 7" size allows easy and immediate display of the menu steps and an intuitive color identification of the alarms.

The controller is equipped with a system of **3 independent microprocessors** that communicate with each other using **CAN-BUS** connection system, which is well know for known for its high reliability and immunity to electrical interference.









## LATEST GENERATION controller

### COLD CONNECTIVITY AND INTELLIGENT CONTROL

The ACP7 controller has innovative features that reduce consumption and thanks to 3 separate microprocessors for: REGULATION, ALARM and DISPLAY allows continuous monitoring of parameters.

The **REAL-TIME** function allows daily display of temperatures. The functionality is completed by the events lists. Temperature recording and operating variables sampled every 30 sec. for up to **10 years.** All data is stored and the event log files can be downloaded to a PC via the built in **USB port.** 

The **ECOMODE** function enables the ULT freezer to operate in energy saving mode.


It is possible to regulate access to the controller and enable the modification of the user parameters through the use of **PERSONAL PASSWORDS**.

The **SMART DIAGNOSTIC** function ensures constant monitoring of the wear of the main components. When the preset maximum wear threshold has been reached, an alarm will prompt the replacement of the component, reducing the risk of machine downtime.

The internal temperature is constantly monitored with an alert when it reaches the safety threshold.

The **BACK-UP BATTERY** provides periodic tests, state of charge and replacement warning. In case of a power failure, the battery has a 30 hours operating time.

In case of power or mechanical failure, the **BACK-UP CO2** or **LN2** system (optional) guarantees the safety of the samples for approximately 24 hours. In case of prolonged door opening, the CO2 or LN2 back up system starts working to support the immediate recovery of the temperature, thanks to the intelligent reading of the electronic system.





AUTO







MORE SAFETY FOR YOUR SAMPLES






# LATEST GENERATION controller

### CONSTANT MONITORING, WHERE AND WHEN YOU WANT IT

The new controller ensures maximum connectivity and constant monitoring thanks to the **CloudALS** computing, or local storage server (Owner Server), easily accessible from PC, Smartphones and Tablets.

Using this option, users can check, in Real-Time, the storage temperatures of the ULT freezer. The service engineers can promptly access detailed information on the operation of the refrigerating unit and its components.

The possibility of choosing between a **3G** network, or **Ethernet**, or **Wi-Fi** communication module facilitates any type of connection in the work environment, allowing the remoting of the functions of one or more Angelantoni devices (refrigerators, ULT freezers, etc.) that has the latest generation controller (ACP series) and that are installed in a single, or in several laboratories.









#### **REMOTE FUNCTIONS**

- Master Data management
- > Multi-Track Real-Time graphics
- > E-mail events notification

JSER

- Real-Time variable consultations
- Reading of configuration parameters
- Writing of configuration parameters
- Firmware updating
- > Teleservice

9

SERVICE





## **OPTIONALS**

OFTIONALS		NEXUS						NEXUSSIim	
	340V-3/340SV-3	340V-4/340SV-4	530V-3/530SV-3	530V-4/530SV-4	110H/110SH	370H/370SH	550H/550SH	810V-4/810SV-4	810V-5/810SV-5
7-days temperature chart recorder for SH and SV models					11351				
7-days temperature chart recorder for H and V models					11352				
BIOGUARD Card door openining system for one door (only when ACP7 istalled)	14161			N.D.	14161				
N.1 Additional user badge for BIOGUARD system					11976				
Porth hole 23 mm with cap					11234				
SS internal shelf (Max 10)	112	04	11205		N.D.			13006	
	112	02	11203		N.D.			13004	
SS sliding telescopic drawer	max9 drawers	max 8 drawers	max 9 drawers	max 8 drawers	N.D.	N.D.	N.D.	max 12 drawers	max 10 drawers
CO2 back-up system (up-right models V)	14163								
CO2 back-up system (chest models H)	14164								
LN2 back-up ystem with hose (all models H e V)	14165								
GSM Module	13884								
Additional PT 100 probe (three wires), class A	13233								
Additional PT 100 class A probe with 4-20A signal converter	11870								
Addit. cert. PT 100 probe,class A with 3 points temp. Calib. (-30°C,-40°C,-43°C) only for model SV e SH	13665								
Addit. certified PT 100 probe, class A with 3 points temp. Calib. (-60 $^{\circ}$ C/-70 $^{\circ}$ C/-80 $^{\circ}$ C) only for model H e V					14174				
Ethernet Module	13921								
CRYOLOG License	OG License 13925								
ALS IQ and OQ protocols					IQ_OQ ALS				

### Rack for cryoboxes 50mm

Model	Code	Туре	Material	Racks Tot.	Boxes Tot.	Vials Tot. (1ml)		
	14057	Drawer	Stainless Steel					
Nexus 340 V-3	14067	Llive	Stamess Steer	9	252	25200		
	14069	Hive	Aluminium					
	14057	Drawer	Stainless Steel					
Nexus 530 V-3	14067	Hive	Stanless Steel	12	336	33600		
	14069	IIIve	Aluminium					
	14014 (1)	Drawer	Stainless Steel	10				
	14463 <mark>(2)</mark>	Diawei	otalilless oteel	2				
Nexus 340 V-4*	14019 (1)		Aluminium	10	230	23000		
Nexus 540 V-4	14465 <mark>(2)</mark>	Hive	Aldminium	2	230			
	14016 (1)	піче	Stainless Steel	10				
	14464 <mark>(2)</mark>			2				
	14014 (1)	Drawer	Stainless Steel	15	315	31500		
	14463 <mark>(2)</mark>	Diawei		1				
Nexus 530 V-4*	14019 (1)		Aluminium	15				
	14465 <mark>(2)</mark>	Hive	Aldining	1				
	14016 (1)		Stainless Steel	15				
	14464 <mark>(2)</mark>		Stanless Steel	1				
	14013	Drawer	Stainless Steel		576			
NexusSlim 810 V-4	14018	Hive	Aluminium	24		57600		
	14015	Inve	Stainless Steel					
	14005	Drawer	Stainless Steel					
NexusSlim 810 V-5	14020	Hive	Aluminium	30	480	48000		
	14017	THVC	Stainless Steel					
Nexus 110H	14497	14497 Alumini		9	72	7200		
	14494	inve	Stainless Steel		12	7200		
Nexus 370H	14010	Hive	Aluminium	21	231	23100		
	14008		Stainless Steel		201	20100		
Nexus 550H	14010	Hive	Aluminium	33	363	36300		
110,43 00011	14008		Stainless Steel		000			



1 - RACK WITH 20 BOXES 2 - RACK WITH 15 BOXES

For these two types of models it is need to order both types of racks, for example if beehive type material is chosen, you must order rack code



Code	CARDBOARD BOXES
14021	Cardboard white standard box (133x133x50), with divider 9x9 for vials 2 ml
14331	Cardboard white standard box (133x133x50), with divider 10x10 for vials 1 ml





### MODELS

HORIZONTAL								
-86°C		-40°C						
Model	Code	Model	Code					
NEXUS 110 H	14082	NEXUS 110 SH	14149					
NEXUS 370 H	14083	NEXUS 370 SH	14150					
NEXUS 550 H	14084	NEXUS 550 SH	14151					
VERTICAL								
-86°C -40°C								
Model	Code	Model	Code					
NEXUS 340 V-3*	14086	NEXUS 340 SV-3*	14153					
NEXUS 340 V-4*	14088	NEXUS 340 SV-4*	14154					
NEXUS 530 V-3*	14087	NEXUS 530 SV-3*	14155					
NEXUS 530 V-4*	14089	NEXUS 530 SV-4*	14156					
NEXUSSIim 810 V-4*	14092	NEXUSSIim 810 SV-4*	14159					
NEXUSSIim 810 V-5*	14093	NEXUSSIim 810 SV-5*	14160					
SDECIFICATIONS *Inner doors quantity								

## **SPECIFICATIONS**

MODELS				HORIZONTAL		VERTICAL				
				NEXUS		NEX	KUS	NEXUSSIim		
			110 370		550	340 530		810		
DIMENSIONS										
Ext. Dim. (mm) (WxDxH) Int. Dim (mm) (WxDxH)			730 x 990 x 1100 450 x 560 x 450	1780 x 920 x 1140 1050 x 520 x 670	2320 x 920 x 1140 1590 x 520 x 670	791 x 1026 x 1978 450 x 530 x 1350	1040 x 1026 x 1978 700 x 530 x 1350	1112 x 1114 x 1978 906 x 601 x 1472		
Liters (It)			113	370	550	322	500	801		
Weight (kg)			150 (SH) 234 (SH) 180 (H) 260 (H)		275 (SH) 312 (H)	214(SV) 262(SV) 234(V) 300(V)		312 (SV) 362 (V)		
TECHINCAL PERFORMANCES										
Range of temperature (°C)	н	V	(- 40 ~ - 86)							
Range of temperature ( C)	SH	sv	(-20 ~ -43)							
Working temperature (°C)	н	V	- 80							
SH SV			- 40							
Noise level (db (A))					1	< 52				
Maximum absorbed current (A)		V	7	7	7	7	9	9		
	SH	SV	3,5	3,5	3,5	3,5	3,5	4		
Power supply					220 V	– 115 V (50 Hz o 60	Hz)			
Refrigerant Gas	V	MIX (R23+R290+R170)(for -86°C)								
	SV	R452a (for -40°C)								
STRUCTURE										
Exterior finish			Powder coated epoxy paint							
Internal material			Stainless steel AISI 304							
Insulation (mm)					140 (PU)			70 (VIP)		
Inner doors (nr)			1	1 4 4			or 4	4 or 5		
Shelves (nr)				1		2 or 3 heigh	t-adjustable	3 or 4 height-adjustable		
Max load per shelf (kg)				1		40 50				







# **ALS** SERVING LIFE SCIENTISTS



**Angelantoni Life Science (ALS)** is a world leader in the supply of refrigeration equipment and in the design of technological solutions in the biomedical sector, constantly engaged in innovation and safety, both biological and environmental.

**AS** brand provides for a wide range of refrigerators cabinets, ultra-low temperature freezers, blood banks refrigerators, freezers to preserve blood components, mortuary prefabricated rooms, stability tests and plant growth chambers, refrigerators and freezers for COVID-19 vaccines.

**AG** brand supplies standardized brine chilling units for pharmaceutical and chemical applications or pharma process, manual or automated solutions allow to reach -70°C for special walk-in chambers and shelters to storage vaccines or other farmaceutical products.

**STERIL** brand provides for equipment able to meet any product protection need, as well as the product, operator and environment safety requirements, for any level of concentration and for any kind of substance (horizontal and vertical laminar flow cabinets, biohazard and cytostatic safety cabinets, laminar flow pass boxes with UV, sanitized hydrogen peroxide pass boxes, weighing, sampling and dispensing cabinets and isolators designed in accordance with the latest international standards (GMP).

**AIC** brand provides Waster, and an automated treatment system for hospital and contaminated waste.





For more information, visit www.cik-solutions.com or contact us at info@cik-solutions.com

Angelantoni TECHNOLOGY FOR LIFE